



SF41-series Pillow Speakers

IL824 Specification Sheet

Tek-CARE®400 P5 Nurse Call System • Rev. 7 - 11/2012

Features

- UL® 1069 Listed
- Oxygen safe
- Fully sealed for easy cleaning
- Flame-retardant (UL® 94V-0) high-impact plastic
- Long lasting switches surpass UL®1069 requirements
- 8-foot cord with DIN plug and security clip
- Provides pillow speaker intercom and light control (when used with the appropriate modules)
- Special Order: RCA or Philips TV controls, analog volume control, 10-foot cord, earphone jack

Description

SF41-series Pillow Speakers are used with Tek-CARE®400 IR411- and IR412-series patient stations, Tek-CARE®400 P5 IR411P5 and IR412P5 patient stations, and Tek-CARE® NC300II IR319- and IR320-series patient stations.

The pillow speakers provide call placement to the master station, as well as remote control for hospital-grade digital TVs (Zenith is standard; RCA and Philips are available by special order). When used in conjunction with the appropriate TekTone® pillow speaker modules, the pillow speakers also supply two-way, hands-free communication with the master station, and automatic silencing of TV audio while the intercom is in use.

Push buttons are provided to call the nurse, plus TV power, channel selection up/down, Closed Captioning, mute and two light controls. Buttons are also provided for volume up/down (analog volume controls are available by special order), and for direct channel access (in some models).

Listings



UL® 1069 listed; UL® File Number E165723.



SF41ZKL (standard)

for Zenith Digital TVs, includes direct channel access, digital volume controls and 2 light controls



SF41ZDL (standard)

for Zenith Digital TVs, includes channel up/down, digital volume controls and 2 light controls

SF41-series Pillow Speakers

Required Components

Hospital-grade Zenith Digital TV

for Tek-CARE®NC300II Nurse Call Systems:

- IR319- or IR320-series Patient Station
- PM311C or PM312C Pillow Speaker Intercom Module
- PM321B or PM322B Light Control Module

for Tek-CARE®400 Nurse Call Systems:

- IR411- or IR412-series Patient Station
- PM421 or PM422 Pillow Speaker & Bed Interface Module

for Tek-CARE®400 P5 Nurse Call Systems:

- IR411P5 or IR412P5 Patient Station
- PM421P5 or PM422P5 Pillow Speaker Module

www.tektone.com

277 Industrial Park Road • Franklin, NC 28734 • tektone@tektone.net

Phone: (828) 524-9967 • Fax: (828) 524-9968 • Sales: Choose option 2 • Tech Support: Choose option 3
TekTone®'s quality system is registered by UL® to the ISO 9001 standard. (Reference #10001510.)